

# Physical Properties

Property	Value	Tolerance
Oxygen Permeability	200	$\pm 40$
Specific Gravity	1.18	$\pm 0.002$
Hardness (Shore D)	79	$\pm 2$
Flexural Modulus, MPa	1194	$\pm 200$
Index of Refraction	1.430	$\pm 0.002$
Water Content (%)	$\leq 1$	$\pm 0.5$
Luminous Transmission (380nm-780nm)	87	$\pm 5$
Wetting Angle (advancing/receding)	98/48	$\pm 10$

## Product Selection

Diameter	Clear with UV	Ice Blue with UV	Blue with UV	Green with UV	Blue non-UV	Green non-UV	Orange non-UV
12.7mm		●	●	●	●	●	●
17.0mm	●	●					
19.0mm	●	●					
21.0mm	●	●					
25.0mm	●						